

# Why DNA Quality Matters – From Research to GMP

Dr. Martin Schleef - Founder of PlasmidFactory GmbH

# Agenda

## **1. DNA as the Foundation of Cell & Gene Therapy**

Role of DNA vectors in modern CGT development

## **2. Plasmid Design Challenges**

How sequence and vector architecture influence performance

## **3. Manufacturing of High-Quality DNA**

Key steps from fermentation to purified plasmid DNA

## **4. Quality Control and Product Robustness**

Ensuring consistency from research to GMP

## **5. Advanced DNA Vector Technologies**

Minicircle DNA and next-generation vector formats

## **6. PlasmidFactory's Expertise**

Supporting CGT programs from early research to clinical development



# DNA as the Foundation of CGT

Role of DNA vectors in modern CGT development



# High Quality DNA by PlasmidFactory

PlasmidFactory supports researchers and therapy developers

- ✓ Robust plasmid production
- ✓ Advanced vector formats
- ✓ Seamless transition from research to GMP

25+ years experience | 3500+ plasmid projects | CGT focus



Company info



# The CGT Toolbox

Cell and gene therapy development relies on several key components

## ✓ Cells

- Patient-derived cells
- Production cell lines

## ✓ Gene transfer technologies

- Electroporation
- Lipid nanoparticles, PEI, ...
- Viral delivery

## ✓ Vectors

- AAV
- Lentivirus
- Plasmid / Minicircle DNA
- RNA



**All vector systems ultimately depend on high-quality DNA.**

# DNA Vector Formats

## Common non-viral DNA formats

### ✓ Plasmid DNA

- Most widely used backbone

### ✓ Linear DNA fragments

- Proprietary MIDGE® vectors

### ✓ Minicircle DNA

- Backbone-free expression vectors

### ✓ Vector design strongly influences:

- Transfection efficiency
- Safety
- Manufacturing feasibility
- Regulatory acceptance



# DNA Design Challenges

How sequence and vector architecture influence performance



# Plasmids Are Not Equal

## ✓ Plasmids widely differ in:

- Size
- Sequence complexity
- Topology
- Regulatory elements

## ✓ These differences strongly affect:


- Scalability
- Stability
- Manufacturability
- Biological performance

### Structure in Textbooks

*Ori* = origin of replication

*AB<sup>R</sup>* = selection marker


*GOI* = gene of interest



# Topology of DNA Vectors

## ✓ Plasmids widely differ in:

- Size
- Sequence complexity
- Topology
- Regulatory elements





Lane 1: Dimer form of the Plasmids  
Lane 2: Dimer form *Bam*HI digested  
Lane 3: Dimer form *Hind*III digested  
Lane 4: Mono-/Dimer form of the Plasmids  
Lane 5: Monomer form *Bam*HI digested  
Lane 6: Monomer form *Hind*III digested  
Lane 7: 1 kb ladder (DNA Size Marker)  
(PlasmidFactory; Art. No. MSM-850-50)

Transformation of a 8393 bp Plasmid

### Control digest

*Hind*III: 5417 bp + 2976 bp  
*Bam*HI: 1525 bp + 6868 bp



**FIG. 2** AFM image of the 7-year-stored plasmid DNA. The AFM image shows ccc plasmid DNA of different degrees of condensation and oc forms (without linear forms, chromosomal DNA aggregates, or other contaminants). The figure shows an AFM image with an edge size of 3  $\mu$ m. The altitude

Walther et. al, 2013


# Impact of Plasmid Size

## ✓ Large plasmids create several challenges:

- Reduced bacterial copy number
- Lower fermentation yield
- More complex impurity profile

## ✓ Results:

- More difficult purification
- Higher production cost
- Potential variability between batches



**Expert process design is therefore essential.**

# Sequence Features That Affect Production

## ✓ Certain sequences reduce plasmid stability:

- Inverted terminal repeats (ITRs)
- Long homopolymers (PolyA)
- Hairpin structures
- Repetitive elements



## ✓ These regions can cause:

- Deletions during replication
- Sequence instability
- Reduced vector performance

### Linearization of custom plasmids & Minicircles



### POLYARESCUE®



### ITRPROTECT® / ITRRESCUE®

PlasmidFactory's proprietary technologies enable amplification of plasmid and Minicircle DNA preserving and repairing the sensitive ITR sequences.

C G C C C G G G C A A G C C C G G G C G T C G G G C G A C C T T T G G T C G C C C G  
G C G G G C C C G T T T C G G G C C C G C A G C C C G C T G G A A A C C A C G C C C C

# CpG Content


## ✓ CpG motifs influence plasmid behavior:

- Innate immune activation
- Transcriptional regulation
- Plasmid stability

## ✓ CpG-optimized or CpG-free plasmids can improve:

- Expression levels
- Tolerability in therapeutic applications

CpG influences  
the performance  
of a plasmid



nature  
biotechnology

CpG-free plasmids confer  
reduced inflammation and  
sustained pulmonary  
gene expression

Stephen C Hyde<sup>1,6,8</sup>, Ian A Pringle<sup>1,6,8</sup>,  
Syahril Abdullah<sup>1,6-8</sup>, Anna E Lawton<sup>1,6,8</sup>,  
Lee A Davies<sup>1,6</sup>, Anusha Varathalingam<sup>1,6</sup>,  
Graciela Nunez-Alonso<sup>1,6</sup>, Anne-Marie Green<sup>1,6</sup>,  
Reto P Bazzani<sup>1,6</sup>, Stephanie G Sumner-Jones<sup>1,6</sup>,  
Mario Chan<sup>2,6</sup>, Hongyu Li<sup>3</sup>, Nelson S Yew<sup>4</sup>,  
Seng H Cheng<sup>4</sup>, A Christopher Boyd<sup>5,6</sup>,  
Jane C Davies<sup>2,6</sup>, Uta Griesenbach<sup>2,6</sup>,  
David J Porteous<sup>5,6</sup>, David N Sheppard<sup>3</sup>,  
Felix M Munkonge<sup>2,6</sup>, Eric W F W Alton<sup>2,6</sup> &  
Deborah R Gill<sup>1,6</sup>


Pulmonary delivery of plasmid DNA (pDNA)/cationic liposome complexes is associated with an acute unmethylated CG dinucleotide (CpG)-mediated inflammatory response and brief duration of transgene expression. We demonstrate that retention of even a single CpG in pDNA is sufficient to elicit an inflammatory response, whereas CpG-free pDNA vectors do not. Using a CpG-free pDNA expression vector, we achieved sustained ( $\geq 56$  d) *in vivo* transgene expression in the absence of lung inflammation.

# Selection Markers

Traditional plasmids contain antibiotic resistance genes

## ✓ Implications:

- Regulatory concerns
- Metabolic burden on host cells
- Increased vector size



**Alternative systems and optimized designs can improve vector performance.**

# Consequences for CGT Development

## ✓ Inconsistent DNA quality can lead to:

- Irreproducible research results
- Difficulties in technology transfer
- Delayed clinical development

## ✓ Worst case:

- Late-stage surprises in clinical development



*GMP production facility, PlasmidFactory GmbH, Bielefeld*

**DNA quality must be addressed early in development.**

# Our Approach

✓ **PlasmidFactory addresses these challenges through:**

- Deep expertise in difficult plasmid sequences
- Optimized amplification strategies
- Advanced analytics

✓ **Our goal:**

- Consistent DNA quality from discovery to GMP



# Manufacturing of High Quality DNA

Key steps from fermentation to purified plasmid DNA



# Quality starts with Process Design

DNA quality is not just final QC

## ✓ It depends on:

- Well-defined raw materials
- Validated equipment
- Controlled processes
- Trained personnel

### Pilot Run

- process parameters
- productivity
- specifications



### E. coli Cell Bank

- K12 safety strain
- transformation
- research cell bank
- master cell bank
- working cell bank



**Consistency of manufacturing is the key determinant of quality.**

# Plasmid Fermentation Strategy



Plasmid production uses fed-batch fermentation

## ✓ Advantages:


- Scalable production
- Controlled growth conditions
- High plasmid yield

## Large Scale Fermentation

- bio-reactor
- monitoring of operating parameters
- reproducible process



Steel fermenter




Industrial batches can reach multi-gram DNA scale.

# Plasmid Fermentation Strategy

Steel or single-use

A question of cleaning validation!




R. Shankar et al. 2020  
(PlasmidFactory, not published)

## Fermentation in Single-Use Fermenters for plasmid and minicircle DNA production

F Weber<sup>1</sup>, C Buschmann<sup>1</sup>, R Shankar<sup>1</sup>, R Baier<sup>1</sup>, C Ramke<sup>2</sup>, R Klementz<sup>2</sup>, M Leupold<sup>2</sup>, J Rupprecht<sup>2</sup>, A Vetter<sup>3</sup>, K Teschner<sup>3</sup>, S Klausning<sup>3</sup>, M Schmeer<sup>1</sup>, M Schleef<sup>1</sup>

1: PlasmidFactory GmbH, Meisenstrasse 96, 33607 Bielefeld, Germany | 2: Sartorius Stedim Biotech GmbH, August-Spindler-Straße 11, 37079 Göttingen, Germany | 3: Sartorius Xell GmbH, Waldweg 21, 33758 Schloss Holte-Stukenbrock, Germany



SARTORIUS

# Cell Lysis and Clarification

Efficient lysis is required to release plasmid DNA

## ✓ Key considerations:


- DNA integrity
- Removal of host cell impurities
- Process scalability

## Cell Harvest and Alkaline Lysis

- bulk biomass
- cleared lysate



## In-Line Lysis to replace bottle shaking



(12) <b>United States Patent</b> Voss et al.	(10) <b>Patent No.:</b> US 7,842,481 B2
	(45) <b>Date of Patent:</b> Nov. 30, 2010
(54) <b>METHOD FOR PRODUCING EXTRA-CHROMOSOMAL NUCLEIC ACID MOLECULES</b>	(56) <b>References Cited</b>
(75) <b>Inventors:</b> Carsten Voss, Lage (DE); Erwin Flaschel, Bielefeld (DE)	U.S. PATENT DOCUMENTS
(73) <b>Assignee:</b> Plasmid Factory GmbH + Co. KG, Bielefeld (DE)	5,407,426 A * 4/1995 Spears ..... 604:24
	2005/0026177 A1 * 2/2005 Urhaker et al. .... 435:6
	2005/0079534 A1 * 4/2005 Warthen et al. .... 435:6
	FOREIGN PATENT DOCUMENTS
	WO WO 03102184 A1 * 12/2003

Modern processes often use in-line lysis systems.

# Purification Strategy & Impurity Removal

## ✓ Purification typically combines:

- Capture chromatography
- Polishing chromatography
- Tangential flow filtration

### Chromatography

- procedure depends on required quality
- removes undesired plasmid forms and chromosomal DNA



### Fill and Finish


- purified plasmid DNA
- aseptic filling




## ✓ Major impurities include:

- Chromosomal DNA
- RNA
- Host cell proteins
- Endotoxin

Agarose gelelectrophoresis



SDS-PAGE





6.1 kb plasmid

- 1: Size standard
- 2: Lysate
- 3: Lysate filtered
- 4: Eluate Capturing chromatography
- 5: Eluate Polishing chromatography
- 6: Eluate Polishing 2 chromatography
- 7: Eluate Polishing EF chromatography
- 8: Desalted by TFF
- 9: Filtered product

## ✓ Optimized purification strategies ensure:

- High purity
- Consistent batch quality



I. Schmitt, R. Shankar, et al., not published

**Goal: High purity plasmid DNA suitable for therapeutic applications.**

# Quality Control & Product Robustness

Ensuring consistency from research to GMP





# Monitoring Plasmid Topology

Plasmid topology strongly influences biological activity

## ✓ Typical forms:

- Supercoiled
- Open circular
- Linear




**Validation of an analytical capillary gel electrophoresis method for the relative quantification of plasmid DNA topoisomers**

A Burgardt, M Strakeljahn, R Shankar, M Wiedemann & M Schleaf

## CGE Specificity

**Validation of an analytical capillary gel electrophoresis method for the relative quantification of plasmid DNA topoisomers**


A Burgardt, M Strakeljahn, R Shankar, M Wiedemann & M Schleaf



## CGE Linearity

**Validation of an analytical capillary gel electrophoresis method for the relative quantification of plasmid DNA topoisomers**


A Burgardt, M Strakeljahn, R Shankar, M Wiedemann & M Schleaf



## CGE Repeatability

**Validation of an analytical capillary gel electrophoresis method for the relative quantification of plasmid DNA topoisomers**

A Burgardt, M Strakeljahn, R Shankar, M Wiedemann & M Schleaf



Analytical methods such as capillary gel electrophoresis (CGE) allow precise monitoring.


# Product & Freeze/Thaw Stability

Plasmid DNA is generally stable under appropriate conditions

## ✓ Examples:

- Storage in WFI at suitable concentration
- Storage at  $-20\text{ }^{\circ}\text{C}$  for long-term preservation

## Long term stability of plasmid DNA




Plasmid DNA in WFI  
at  $1,0\text{ mg mL}^{-1}$   
Storage at  $-20\text{ }^{\circ}\text{C}$

Repeated freeze-thaw cycles have limited impact on plasmid integrity.

## Impact of freeze-thaw-cycles on plasmid DNA $\geq 10\text{ kb}$

Number of cycles over 2 months



Number of cycles	ccc peak [%]	oc peak [%]
2	98,8	1,1
8	98,5	1,6
16	99,1	0,9

## ✓ This robustness is important for:

- Distribution
- Storage
- Routine laboratory use
- Clinical applications

Long-term stability supports reliable use in research and clinical applications.

# Advanced DNA Vector Technologies

Minicircle DNA and next-generation vector formats



# Beyond Classical Plasmids

To address limitations of conventional plasmids, alternative DNA formats are used:

## Main vector concepts:

- ✓ **Plasmid DNA (circular, with backbone):**
  - Standard format for gene delivery
  - Flexible, proven, but includes bacterial elements
- ✓ **Backbone-free circular DNA (Minicircle):**
  - Circular expression cassette only
  - No bacterial sequences, optimized for expression



## Additional approaches:

- ✓ **Synthetic and enzymatic linear DNA (e.g. dbDNA, MIDGE<sup>®</sup>):**
  - Backbone-free and cell-free production possible
  - Covalently closed linear formats improve stability
  - Limitations in scalability, handling and expression performance depending on format
- ✓ **Other "mini" or "nano" plasmids:**
  - Reduced bacterial backbone
  - Partial improvement, but not fully backbone-free

**Alternative formats to plasmids can improve safety and efficiency.**

# Linear DNA: MIDGE® Vectors

MIDGE (Minimalistic Immunogenically Defined Gene Expression) vectors are linear DNA constructs without bacterial backbone elements.

## ✓ Key characteristics:

- Linear DNA expression cassette
- Absence of bacterial sequences
- Reduced vector size compared to plasmids


## ✓ Potential advantages include:

- Lower immunogenicity
- Improved safety profile
- Simplified vector design

## ✓ Potential disadvantages (incl. dbDNA):

- No enzyme-free
- Not supercoiled

Linear DNA (MIDGE) produced *in vitro* by use of enzymes?




MIDGE - a registered trademark of PlasmidFactory

# Minicircle: Advantages at a Glance

## Bacterial backbone-free Minicircles

- ✓ No antibiotic resistance or other bacterial markers
- ✓ Less immunogenicity
- ✓ Lower DNA toxicity
- ✓ Higher transfection efficiency
- ✓ Reduced transgene silencing
- ✓ Stronger, more stable gene expression
- ✓ Almost no cargo size restriction
- ✓ Improved yield – minimized costs
- ✓ No backpacking in AAV production
- ✓ Ideal for virus-free gene transfers

**Successful in several clinical trials**



Kobelt et al, 2012




Minicircle info




# Minicircle Technology: Proprietary Production Method

Minimalistic DNA vector reduced to the gene of interest (GoI)

## Minicircle production



## Minicircle vs. plasmid features



# Minicircle Purification & Quality Control

## ✓ Cultivation:

- Recombinase expression
- Unidirectional
- Extremely fast


## ✓ Purification:

- Separation of Minicircle and parental plasmid
- Chromatography-based process (e.g. affinity)
- Removal of backbone and process-related impurities


## ✓ Quality control:

- Analysis by capillary gel electrophoresis (CGE)
- Verification of monomeric supercoiled topology
- Confirmation of purity and residual MP removal


**Arabinose induction of recombination:**  
Extremely fast process with high efficiency





**Separating MP and MC:**  
Affinity chromatography



**Product Quality determined by CGE:**  
monomeric, supercoiled, removal of chromosomal DNA



# Clinical Scale Generation of functional CAR-T Cells using a Minicircle-based Sleeping Beauty Transposon System



**Díez et al., 2025, Mol Ther Methods Clin Dev**



- ✓ MC approach results in equal tumor eradication to LVV (A)
- ✓ MC shows improved genomic safety vs. LVV (B)
- ✓ MC approach results in prolonged survival equal to LVV (data not shown)
- ✓ MC un-detectable in final CAR T product (data not shown)

 [Article title:](#) Generation and GMP scale-up of human CAR-T cells using non-viral Sleeping Beauty transposons for B cell malignancies

GMCT Webinar



# Minicircles in AAV Production for higher Vector Purity and increased Titers



**Kraemer et al., ESGCT Annual Congress 2021**  
**Schnödt et al., Mol. Ther.-Nucl. Acids 2016**

- ✓ **Marked increase in viral titers:** 3 to 4-fold (A)
- ✓ **No back-packaging:** higher vector purity and product consistency (B)
- ✓ **Improved transduction efficiency:** up to 30× higher (scAAV; data not shown)
- ✓ **Cleaner & safer** – lower immunogenicity, less silencing, up to GMP Grade (data not shown)


 [Article title: DNA Minicircle Technology Improves Purity of Adeno-associated Viral Vector Preparations](#)


AAV webinar



# Keeping ITRs Intact: ITRPROTECT® & ITRRESCUE®

Preserved ITRs for better yield, purity and packaging efficiency



- ✔ **ITR instability is a known challenge:** In AAV transfer plasmids, inverted terminal repeats (ITRs) are prone to truncation during plasmid amplification, reducing vector yield, purity, and packaging efficiency ([Radukic et al., Mol. Ther.-Nucl. Acids 2025](#)). 
- ✔ **ITRPROTECT® safeguards your sequences:** Our proprietary technology ensures stable amplification while preserving ITRs exactly as provided by the customer.
- ✔ **ITRRESCUE® restores functionality:** We can analyze ITR integrity and, if necessary, replace truncated regions with full-length wild-type ITRs to recover fully functional AAV plasmids.
- ✔ **Reduce variability:** Full-length ITRs, either restored when needed and preserved during amplification, ensure stable vector performance and predictable AAV packaging.

# PlasmidFactory's Expertise

Supporting CGT programs from early research to clinical development



# Your trusted CDMO Partner



## Why chose PlasmidFactory



- ✓ 25+ years of DNA manufacturing expertise
- ✓ Independent, flexible & trusted European CDMO
- ✓ Scalable production: Research → High Quality → GMP
- ✓ Experts in cell & gene therapy and vaccination apps
- ✓ Experts for difficult sequences (ITRs, large plasmids, poly(A), repetitive & CpG-rich regions)
- ✓ Proprietary technologies: Minicircle, ITRPROTECT<sup>®</sup>, ITRRESCUE<sup>®</sup>, POLYARESCUE<sup>®</sup>, MIDGE<sup>®</sup>
- ✓ 3,500+ plasmid & Minicircle projects delivered
- ✓ 99.9% success rate, proven in clinical trials



# Available Quality Grades




## Scientific Quality (SQ) Grade

Advanced research standard

- ✓ 4 parallel production lines
- ✓ Scalable process starting from 0.5 mg\*
- ✓ 2 options available:
  - *Research Grade* for basic requirements
  - *CCC Grade* with ≥ 95% supercoiled DNA

### Key features

- Bioreactor fermentation

 2 x 30L  
2 x 20 L


## High Quality (HQ) Grade

Highest non-GMP standard

- ✓ 2 HQ facilities
- ✓ According to EMA guidelines
- ✓ Production scale 10 mg – 10 g\*\*
- ✓ Also available as "HQ in GMP" Grade

### Key feature

- Complete traceability

 2 x 30L  
2 x 200L


## GMP Grade

Clinical & commercial compliance

- ✓ Dedicated state-of-the-art GMP facility
- ✓ GMP / FDA-compliant
- ✓ End-to-end single-use processes
- ✓ Late-stage clinical trials and commercial supply

### Key feature

- Single use equipment

 2 x 40L

\*0.5 mg for MC; 5 mg for plasmid DNA

\*\*1 g for MC; 10 g for plasmid DNA; more on request

# GMP Manufacturing You Can Trust



Designed for highest plasmid & Minicircle DNA quality requirements



## State-of-the-art GMP facility

- ✓ Dedicated building for GMP Grade plasmid & Minicircle DNA
- ✓ 450 m<sup>2</sup> cleanrooms suites (grades C & D)
- ✓ Facility designed to prevent cross-contamination

## Flexible manufacturing with highest safety standards

- ✓ End-to-end single-use USP- and DSP- process
- ✓ Ph. Eur. / USP grade WFI supply
- ✓ EU GMP (Annex 1) / FDA compliant aseptic fill & finish
- ✓ QP batch release in-house

## Full GMP compliance

- ✓ Data integrity measures and systems in place
- ✓ Annex 11 / 21 CFR Part 11
- ✓ EU GMP-Part 2 and AMWHV

GMP info



# Minicircle DNA: Clinical Footprint



## Snapshot of registered and planned clinical trials using Minicircle for cell engineering

### Registered clinical trials

2020

#### CARAMBA-1 - SLAMF7 CAR-T (Multiple Myeloma)

Recruitment ended

- Phase I/IIa (FiH) • EUCT 2024-512643-23-00 • [NCT04499339](#)
- Virus-free Sleeping Beauty gene transfer using Minicircle DNA + transposase mRNA
- Opened Jun 2020 • 33 participants • DE site (multi-country approvals reported in EU project updates)

2024

#### TranspoCART19 - CD19 CAR-T (r/r B-cell Lymphoma)

Recruiting

- Phase I/IIa • EUCT 2024-514544-90-00 • [NCT06378190](#)
- Non-viral Sleeping Beauty transposon delivered as a Minicircle + SB100X mRNA
- Start Mar 2024 • 27 participants • 9 sites (Spain)

2025

#### LION-1 - ROR1 CAR-T (ROR1<sup>+</sup> tumors)

Recruiting

- Phase I (FiH) • EUCT [2024-512019-36-00](#) • NCT pending (CTIS)
- Start Apr 2025 • 46 participants • 1 site (Germany)
- Includes hematologic + solid tumors (e.g., MCL, CLL, TNBC, OC, ACC)

### Studies in preparation

#### Oncology (engineered cells)

- ROR2 CAR-T (Uniklinikum Würzburg)
- FLT3 CAR-T (Uniklinikum Würzburg)
- SLAM/F7 program (T-CurX)

#### Autoimmune / gene therapy

- Nimvec™ AM510 (Amarna; Type 1 Diabetes)

#### Undisclosed / confidential

- At least 4 additional clinical studies in preparation

# Summary & Take-Home Message

✓ **Successful cell & gene therapy development requires:**

- High-quality vector design
- Robust manufacturing processes
- Validated analytics & QC

✓ **PlasmidFactory supports this with:**

- Expertise
- Innovation
- Full GMP grade manufacturing



SUMMARY



# Thank you for your interest!

Feel free to contact:

[Martin.Schleef@plasmidfactory.com](mailto:Martin.Schleef@plasmidfactory.com)

## PlasmidFactory GmbH

Meisenstraße 96 | 33607 Bielefeld | Germany

[www.PlasmidFactory.com](http://www.PlasmidFactory.com)

